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
General remarks about the practical tests

DEAR PRATICIPANS 參賽者

The practical test are organized in four different laboratories. 實作在 4 間不同實驗室進行

Nº 1- Plant Anatomy, Systematics and Physiology 植物解剖學，分類學和生理學

Nº 2- Animal Anatomy, Systematics and Ecology 動物解剖學，生態學和分類學

Nº 3- Biocheminstry 生物化學

Nº 4- Microbiology 微生物學

- You have **1 hour** in laboratories Nº 1 and Nº 2. 實驗 1 及 2 為 1 個小時
- You have **1 hour 30 minutes** in laboratories Nº 3 and Nº 4. 實驗 3 及 4 為 1.5 個小時
- You can score maximum **40 points** in each laboratory, which means a total of **160 points** for the practical test. 每個實驗你最多可得 40 分，實作試驗共 160 分

Good luck !!!!!!!

Practical test N° 3: Biochemistry

Enzymatic determination of glucose

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TASK 1: You have to perform a calibration curve using a standard of glucose, with known concentration. Then, plot the results as absorbance versus glucose concentration (15 points)

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Important: Raise the red card to call the lab assistant when you are ready to use the spectrophotometer.

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Introduction: 內容介紹

Glucose oxidase (GOD) catalyzes the oxidation of (beta)-D-glucose to D-gluconic acid and hydrogen peroxide. It is highly specific for (beta)-D-glucose and does not act on (alpha)-D-glucose. The horseradish peroxidase (POD) breaks down hydrogen peroxide into water and oxygen, using the dye (4-aminophenazone) as an electron donor. At the same time, the dye is converted to its oxidized form, which is a colored compound. Since the amount of hydrogen peroxide produced indicates how much reaction has taken place, the formation of the red color can be used to follow the course of the reaction.

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
Its major use is in the determination of free glucose in body fluids. Although specific for (beta)-D-glucose, glucose oxidase can be used to measure total glucose, because as a result of the consumption of (beta)-glucose, (alpha)-glucose from the equilibrium is converted to the (beta)-form by mutarotation.

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PRINCIPLE 實驗原理

The reaction system is as follows: 其反應式如下圖

Glucose oxidase reagent: solution containing glucose oxidase, peroxidase, 4-aminophenazone (4-AP), and phosphate buffer pH 7.0 containing hydroxybenzoate.

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
Counter (容量計)

Adjustment method 調整方法

You have to pull up the Adjustable Wheel, then you can revolve the adjustable wheel or knob. Adjust the required volume and push down the Adjustable Wheel.

Remember that minimal and maximal volumes for P100 are 10 μ l and 100 μ l respectively.

For P1000 minimal volume is 100 μ l and maximal volume is 1000 μ l.

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Usage method: 使用方法

Please secure the suction tip, after that slightly push down the pushing knob to first stop, hold and immerse the tip into solution vertically. The immersed depth of the tip is 2-4 mm, then release the pushing knob slowly and make it return to the original position. Take off the pipette from the liquid and place the suction tip of the pipette into a special container receiving the dispensed liquid. The tip must be close to the inner wall of the container. Depress the pushing knob to the first stop and further more to discharge the solution completely from the tip. After that, you can take away the pipette and release the button. Eject the used tip to the trash recipient by pressing the Tip ejecting knob.

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EXPERIMENTAL PROCEDURE 實驗流程

1) Label five 1.5-ml microtubes 1/2 through 1/32 with a marker pen. Using the glucose standard solution ($5 \text{ mg} \cdot \text{ml}^{-1}$) perform the following serial dilutions (in distilled water) in a final volume of 100 μl : 1/2, 1/4, 1/8, 1/16, and 1/32.

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2) Mix well and perform (in a new 1.5-ml microtubes set) the enzymatic determination of glucose for each dilution according to the following scheme.

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3) Mix well and incubate microtubes at 37°C for 5 min.

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4) Put the content of each microtube in a spectrophotometric cuvette.

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5) Read absorbance in a spectrophotometer at 505 nm. **You have to use the blank for calibration. (When you are ready, to read in the spectrophotometer, please call the lab assistant)**

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6) Plot the absorbance versus the amount of glucose (in μ g) on the plotting paper below.

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TASK 2: Determination of the glucose concentration in a sample employing the standard curve obtained before. (10 points)

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1) Perform the glucose oxidase reaction, to the glucose sample (unknown concentration), according to the follow scheme.

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3) Put the content of each microtube in a spectrophotometric cuvette

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4) Read absorbance in the spectrophotometer at 505 nm. You have to use **the blank for calibration. (When you are ready to read in the spectrophotometer please call the lab assistant)**

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5) Using the calibration curve, calculate the glucosa concentration of the sample in $\mu\text{g} \cdot \text{ml}^{-1}$.

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QUESTION 1: As many glucose assays measure the peroxide produced by the glucose oxidase reaction, it is important that the enzyme used for these assays presents: (1.5 points) Deleted 已刪除，不要做。

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WRITE DOWN THE LETTER CORRESPONDING TO CORRECT ANSWER

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Answer:.....Deleted..... 已刪除，不要做。.....

QUESTION 2: Glucose oxidase reagent may contain catalase. If such a condition is not taken into account the obtained results will give **(1.5 points):**

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Answer:.....

QUESTION 3: The most favorable pH value (the point at which the enzyme is most active) is known as the optimum pH. Extremely high or low pH values usually result in a complete loss of enzyme activity due to **(1 point):**

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SELECT ONLY ONE CORRECT ANSWER MARK IT WITH A CROSS.

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- A. B. C. A, B
- A, C B, C. A, B, C.

QUESTION 4: Glucose oxidase from the fungus *Aspergillus niger* was overexpressed in yeast. The glucose oxidase was purified and glycosylation pattern was analyzed by treatment with endoglycosidase H and α -mannosidase. After treatment, an aliquot was used for SDS-PAGE (electrophoresis in polyacrylamide gels containing sodium dodecyl sulphate) in reducing conditions. The remaining enzyme was employed for determination of the K_M (Michaelis-Menten constant) with glucose as the substrate. Michaelis-Menten constant (K_M) is the concentration in moles/litre of a substrate at half the maximum velocity of an enzymatic reaction. **(7 points)**

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The values of the K_M for each glycoforms are shown below the figure 1.

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
Figure 1: Analysis of the deglycosylation of Glucose

Oxidase by 7.5% acrylamide SDS-PAGE gel electrophoresis. Lane 0 is the molecular mass standard. Lane 1 untreated enzyme. Lane 2 endoglycosidase H treated enzyme. Lane 3 α -mannosidase treated enzyme. Lane 4 endoglycosidase H and α -mannosidase treated enzyme (fully deglycosylated enzyme).

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Which of the following is/are correct conclusion(s) from these results obtained in the determination of KM for each glycoform. (figure 1)

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A) Glucose oxidase is a homodimer with a molecular mass of 96 kDa.

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B) The deglycosylated form has a molecular mass of approximately 68 kDa.

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C) Glucose oxidase is glycosylated since the treatment with endoglycosidase H and α -mannosidase results in a form with lower molecular mass.

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D) The polysaccharide moiety of glucose oxidase contains *N*-acetylglucosamine and mannose.

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A B C D

From the results obtained in the determination of K_M for each glycoform, the following conclusion could be made. Which of the following statements is/are correct?

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A) The affinity of the fully glycosylated enzyme for glucose is higher than the affinity of the deglycosylated enzyme

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B) The glucose oxidase activity is completely abolished in the deglycosylated form

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C) The lack of the sugar moiety could cause changes in the structure of the active site of the enzyme resulting in the observed modifications of the K_M s.

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A B C

QUESTION 5: The purified enzyme was analyzed by SDS-PAGE in reducing (DDT+) as well as in non-reducing (DDT-) conditions. The obtained results are shown in figure 2 (4 points)

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
Figure 2: SDS-PAGE analysis of purified glucose oxidase.

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Taking in account the results obtained from figure 1 and figure 2 the most probable conformation of the glucose oxidase is:

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A) A monomeric enzyme non-glycosylated.

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B) A monomeric enzyme glycosylated.

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C) A homodimer consisting of two monomers both glycosylated.

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D) A heterodimer consisting of two subunits one of them glycosylated.

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