Student Code: _____

20th INTERNATIONAL BIOLOGY OLYMPIAD

12th – 19th July, 2009

Tsukuba, JAPAN

PRACTICAL TEST 2

BIOCHEMISTRY

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Total Points: 100

Duration: 90 minutes

Dear Participants,

- In this test, you have been given the following 2 tasks:
 - Task 1: Measurement of acid phosphatase activity (70 points)
 - Task 2: Protein determination (30 points)
- 本實作包括兩部分：
- 第 1 部分：酸性磷酸酶活性的測定 (70 分)
 - 第 2 部分：蛋白質的測定 (30 分)
- **You must write down your results and answers in the ANSWER SHEET. Answers written in the Question Paper will not be evaluated.**
 - Please make sure that you have received all the materials and equipment listed for each task. If any of these items are missing, please raise your hand.
 - At the end of the test, put the Answer Sheet and Question Paper in the envelope. The supervisor will collect this envelope.
 - This series of practices are time consuming. You will need to be well organized and work quickly to complete the five tasks.
 - 你必須把結果及答案寫在答案紙上。寫在考卷上的答案不給分。
 - 確定材料與器材是否齊全，若有缺失，請舉手告知監試人員。
 - 本實作結束後，請將答案紙及考卷一起放入信封中，待監試人員收回。
 - 本實作相當費時，你的操作必須很有組織而且快速，才能完成五個部分。

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Fig. 5

Introduction

Acid phosphatase liberates phosphate ions from phosphorylated molecules under acidic conditions. The purpose of this experiment is to determine the specific activity of the acid phosphatase. You will measure activities of the acid phosphatase using a crude extract from potato in Task 1, and determine a protein concentration of the crude extract in Task 2. Specific activity, which is the activity per unit time per unit weight of protein, is obtained from Tasks 1 and 2. Specific activity is an index of purity; it increases as the enzyme is purified.

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Caution 注意

1. You will be handling small amounts of toxic substances (*p*-nitrophenol and NaOH). You can choose to wear disposable gloves and safety goggles in the experiments if you like.
2. In calculations where answers to previous questions are needed, partials marks will be given if calculated formulas are correct, even if answers are incorrect.
1. 你將會使用少量的有毒物質 (*p*- nitrophenol and NaOH)，如有需要，你可穿戴丟棄式手套及護目鏡。
2. 承上述的實驗結果來進行計算時，只要你計算公式正確，即使答案錯誤，仍可有部分得分。

<u>Materials and Equipments 材料與器材</u>	Quantity
1. Spectrophotometer 光電比色計	1
2. Micropipettes (P1000) 微量吸管	2
3. Micropipettes (P200) 微量吸管	1
4. Tips (one box each for P1000 and P200) 微量 Tips	2
5. Plastic cuvette 塑膠光電比色管	1
6. Test tube holder that accommodates 6-1 to 6-6 可容納六根試管(6-1 to 6-6) 的試管架	1
6-1. Crude extract of acid phosphatase (4 ml in a 15-ml plastic tube, labeled “1X” enzyme) 酸性磷酸酶的粗萃液 (4 ml , 置於 15-ml 的塑膠試管中並標示 1 倍“1X”的酵 素)	1
6-2. 0.5 M Na acetate buffer (pH 5.6) (2 ml in a 15-ml plastic tube) 0.5 M 醋酸鈉緩衝溶液(pH 5.6) (2 ml , 置於 15-ml 的塑膠試管中)	1
6-3. 5 mM pNPP (4 ml in a 15-ml plastic tube)	1
6-4. 0.5 M NaOH (8 ml in a 15-ml plastic tube)	1
6-5. 3% NaCl (10 ml in a 15-ml plastic tube)	1
6-6. Test tubes (Glass)	6

Task 1 (70 points)


Measurement of acid phosphatase activity

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The activity of acid phosphatase is measured by an enzymatic reaction that converts *p*-nitrophenyl phosphate (pNPP) to *p*-nitrophenol (pNP), liberating phosphate. The product, pNP, absorbs light whose wavelength is 400 nm with an absorption coefficient* ($\epsilon_{400 \text{ nm}}$) of $19000 \text{ M}^{-1} \text{ cm}^{-1}$ at extremely alkaline pH. Reaction mixture for an acid phosphatase is slightly acidic. Thus, it must be alkalized for quantification of pNP. In Task 1, you will measure a time course of the reaction and obtain absorbance change per minute that is caused by 1 ml of crude extract. The absorbance change is converted to concentration change by using $\epsilon_{400 \text{ nm}}$. Then, you will calculate a mol number of pNP molecules produced during the reaction by multiplying the concentration change by a volume of sample that is subjected to the measurement of absorbance.

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*What is an absorption coefficient? 什麼是吸光係數?

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I_0 , intensity of incident light 入射光之強度

I , intensity of transmission light 透射光之強度

Absorbance (A) is a physico-chemical property of solution that expresses to what extent a solute absorbs light at a specific wavelength. Absorbance is in proportion to concentration (C) and light path length (L). The constant in the equation is a value characteristic to the solute, and is termed the absorption coefficient (ϵ). Thus, the relationship is formulated as $A = \epsilon C (M = mol\ litre^{-1}) L (cm)$. Absorbance can be converted to concentration, since ϵ is given and L is 1 cm in this experiment. The dimension of ϵ is $M^{-1} cm^{-1}$, because absorbance is an absolute number without units.

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Two enzyme concentrations are to be examined in Task 1. Find the test tube on which '1X enzyme' is labeled, which contains a crude extract of acid phosphatase. Next, find the 15-ml tube that contains 3% NaCl and remove 1 ml of the solution so that the tube now contains 9 ml of 3% NaCl. Add 1 ml of the 1X enzyme solution to it by using a micropipette, which makes '0.1X enzyme' solution. Label the tube as '0.1X'. Next, find 6 empty test tubes. Label each tube with an enzyme concentration and a reaction time as follows.

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0.1X, 1 min

1X, 1 min

Q.1.1. (10 points) First, make an experimental schedule in order to perform all reactions, by describing start (○) and stop (●) signs for each reaction in the table in the Answer Sheet, allowing at least 1 min between the beginning of each reaction. An example for the reaction of '0.1x, 20 min' has been described in the table in the Answer Sheet.

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Q.1.2. (15+10 points) Perform the enzymatic reactions according to the protocol described below and the schedule you made in Q.1.1. Use a new pipette tip in every manipulation. Agitate a mixture by tapping a test tube immediately after an addition. After you perform all the reactions, measure A_{400} of the samples. Write the obtained values in the table in the Answer Sheet, and plot them in the graph. Please note that since water has been used as a blank, the line will not pass through 0 (zero) on Y-axis (origin).

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Protocol for measurement of acid phosphatase activity

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- 1) Mix 0.12 ml of 0.5 M Na acetate buffer (pH 5.6) and 0.24 ml of 5 mM pNPP in a test tube. Start the reaction by adding 0.24 ml of an enzyme solution.
- 2) After the reaction times of 1, 10, and 20 min, respectively, stop the reaction by adding 0.6 ml of 0.5 M NaOH. NaOH stops the reaction and converts the pNP produced into a yellow-colored (A_{400} -absorbing) form.
- 3) After all reactions are stopped, measure A_{400} of the samples.

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Q.1.3. (15 points) Which enzyme concentration gave better linearity in the relationship between time and A_{400} ? Circle the correct one on the Answer Sheet. Read the slope of this straight line from the graph.

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Q.1.4. (5 points) Using the slope obtained in Q. 1.3, calculate the activity in the form of A_{400} change per min per 1 ml of an enzyme solution of concentration “1X”. The length of the light path (L) is 1cm. Your answer should be written with your calculations and the appropriate unit in the Answer Sheet.

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Q.1.5. (5 points) Convert the absorbance change obtained in Q.1.4 to a concentration change by assuming the ϵ_{400} of pNP (to be $19000 \text{ M}^{-1} \text{ cm}^{-1}$). Your answer should be written with your calculations and the unit per min per 1 ml of 1X enzyme solution in the Answer Sheet.

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Q.1.6. (5 points) Convert the concentration change obtained in Q.1.5. to a change in number of moles of pNP. Your answer should be written with your calculations in moles per min per ml of 1X enzyme solution in the Answer Sheet.

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Q.1.7. (5 points) Calculate the total activity (in mol per min) in 4 ml of 1X enzyme solution that was initially given.

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Task 2 (30 points)

Protein determination

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Protein concentration is determined by using a standard protein such as bovine serum albumin (BSA). In Task 2, you will determine a BSA-equivalent concentration of the 1x enzyme solution by the Bradford method. The Bradford method takes advantage of an increase in absorption of Coomassie Brilliant Blue at 595 nm when it is bound to protein.

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By diluting a concentrated BSA solution ($0.4 \text{ mg protein ml}^{-1}$) with 3% NaCl, a 1/2-dilution series was made (0.4, 0.2, 0.1, and $0.05 \text{ mg protein ml}^{-1}$). The dilution series of BSA and the 0.1X enzyme solution, which was made in Task 1, were all similarly treated with dye. Optical density at 595 nm (OD_{595}) of these samples was measured and recorded in the table below.

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Optical density (OD), a measure of the extent to which a substance transmits light or the ‘absorbance’ of suspension of particles.

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Q.2.1.(10 points) Plot OD_{595} against BSA concentration in the graph in the Answer Sheet and depict an approximate straight line.

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Q.2.2.(10 points) Estimate a protein concentration of the 0.1X enzyme solution from the graph, and obtain the protein concentration of the 1X enzyme solution.

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Q.2.3.(10 points) Calculate the specific activity (activity per min per mg protein) of the 1x enzyme solution. Your answer should be written with your calculations and the unit per min per mg protein in the Answer Sheet.

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