Student Code: \_\_\_\_\_

## **22<sup>nd</sup> INTERNATIONAL BIOLOGY OLYMPIAD**

**July 10-17, 2011**

**Taipei, Taiwan**

### **PRACTICAL TEST 1**

### **BIOCHEMISTRY AND CELL BIOLOGY**

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**Total Points: 100**

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**Duration: 90 minutes**

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Dear Participants, 敬愛的參賽者

- In this test, you have been given the following 3 tasks: 本試卷包含三個實做題目  
**Task I:** Protein electrophoresis (35 points) 蛋白質電泳 (35 分)  
**Task II:** Protein quantification (30 points) 蛋白質定量 (30 分)  
**Task III:** Protein purification (35 points) 蛋白質純化 (35 分)
- Check your **Student Code** on the **Answer Sheet** before starting the test.  
檢查答案卷的號碼是否與你的考生編號相同。
- Write down your results and answers in the **Answer Sheet**. **Answers written in the Question Paper will not be evaluated.**  
一定要將答案填入答案卷上的正確位置，答案寫在試卷上將不予計分。
- Make sure that you have received all the materials listed for each task. If any of the listed items is missing, **raise your sign**.  
確認每題的材料與問題是相符的，如果有缺少任何物件，**請舉牌**。
- Use pen only. 不可使用鉛筆
- You should organize your work efficiently but ensure that you complete task II early enough to obtain the spectrophotometer readings to answer the questions that follow.  
你可以根據實際狀況有效率的回答問題，以便完成 task II 中分光光度計的實驗，並回答之後的題目。
- Stop answering **immediately** after the end bell rings. 鈴響後立刻停止作答
- After test, enclose both the **Answer sheets, Question paper, and Data printout** in the provided envelope. Our lab assistants will collect it promptly.  
考試結束後，將**答案卷，試題卷與列印出的數據**放入指定的信封內。助教會迅速的收走。
- No paper or materials should be taken out of the laboratory.  
任何紙張與材料嚴禁攜出實驗室。

Good Luck!!

## Shared instruments: 共用設備

Camera, spectrophotometer, printer 照相機，分光光度計，印表機

## Equipments and Materials: 設備與材料

<u>Equipments: 設備</u>	<u>Quantity</u>
1 Power supply 電源供應器	1
2 Electrophoresis tank (with gel and buffer) 電泳槽 (含膠體與緩衝液)	1
3 Micropipettes P20 and P200 微量吸管 P20, P200	1 each
4 80-well microcentrifuge tube rack 80 孔 微量離心管管架	1
5 Wire test tube rack with 15-mL centrifuge tubes (×6) (yellow cap) 不鏽鋼試管架包括 6 支黃蓋 15 mL 離心管	1
6 4-way test tube rack 四面試管架	1
7 Plastic droppers in 15-mL centrifuge tubes 15 mL 離心管附塑膠滴管	2
8 Micropipette tips (for P20 and P200) 微量吸管頭	1 each
9 Timer 計時器	1
10 96-well microplate 96 孔微量盤	1
11 Marker pen & paper label 奇異筆與貼紙	1 each
12 600-mL beaker for waste disposal 600 mL 燒杯 (裝廢液)	1
13 Scissors 剪刀	1
14 Double-sticker to attach the results 雙面膠 (貼結果用)	1
15 Student Code sticker 學生編號貼紙	1
16 Tissue paper 衛生紙	1
17 Mini centrifuge (if you need to spin down the samples in the microcentrifuge tubes) 微量離心機 (沉澱樣本用)	1

<b><u>Materials: 材料</u></b>	<b><u>Quantity</u></b>
1 Loading dye (microcentrifuge tube-L) (pink tube with orange label) 指示染料 (微量試管-L)(粉紅色試管橙色標籤)	1
2 Pre-stained protein molecular weight marker (microcentrifuge tube-M) (pink tube with orange label) 預染蛋白質分子量指示標記 (微量試管-M)(粉紅色試管橙色標籤)	1
3 Unknown pre-stained protein samples (microcentrifuge tubes-U1 and U2) (pink tube with orange label) 未知預染蛋白質樣本 (微量試管-U1 與 U2)(粉紅色試管橙色標籤)	1
4 CBG reagent in 50-mL centrifuge tube    CBG 試劑裝在 50mL 離心管中	1
5 Bovine serum albumin (BSA) concentration standard (0.5 mg/mL) in microcentrifuge tube (green tube with yellow label) 牛血清白蛋白 (BSA) 濃度測定標準樣本 (0.5 mg/ml) (綠色試管黃色標籤)	1
6 Enzyme E in two microcentrifuge tubes: concentrations X and Y (green tube with yellow label) 濃度 X 與 Y 的酵素 E (綠色試管黃色標籤)	1
7 Distilled water (microcentrifuge tube-ddH <sub>2</sub> O) (green tube with yellow label) 蒸餾水 (微量試管 dd H <sub>2</sub> O) (綠色試管黃色標籤)	1
8 Protein sample (microcentrifuge tube-C) (blue tube with blue label) 蛋白質樣本 (微量試管 C)(藍色試管藍色標籤)	1
9 Anion exchange chromatography column on 15-mL centrifuge tube 15mL 離心管附陽離子交換色層分析管	1
10 Anionic buffers A and B (5 mL each in two separated 15-mL centrifuge tubes) (green cap) 陽離子緩衝溶液 A 與 B (各 5 mL 分別裝於 15 mL 綠色蓋子離心管)	1
11 Coomassie brilliant blue G-250 (CBG) reagent 1 mL in each of six 15-mL centrifuge tubes (A1 to A3 & B1 to B3, red cap) Coomassie brilliant blue G-250 (CBG) 試劑, 分別裝於 15 mL 紅色蓋子離心管, 每管 1mL 並分別標記為 A1, A2, A3, B1, B2, B3)	1

## Task I (35 points)

### Protein electrophoresis 蛋白質電泳

#### Introduction: 簡介

Polyacrylamide gel electrophoresis (PAGE) is a common technique for protein study. It can be used to separate different proteins based on their charges or sizes. A type of PAGE is termed SDS-PAGE, in which the negatively charged chemical, SDS, is added before protein electrophoresis. The amount of SDS that binds to proteins is proportional to the size of the protein which confers each protein a similar charge-to-mass ratio and renders the intrinsic charge of the protein insignificant, at least for this experiment. Thus, the major factor that affects the migration of protein is the molecular weight (MW) of the protein during SDS-PAGE. The relative mobility ( $R_f$ ) of the protein can be calculated as the ratio of the distance migrated by the protein to that migrated by the dye-front. The value of  $R_f$  is inversely proportional to the log of its molecular weight.

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**In the problem set, you will perform the following experiment:**

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1. An electrophoresis tank has been set up for SDS-PAGE, in which a polyacrylamide gel has been secured on electrode assembly and electrophoresis buffer has been filled. There are 10 wells for sample loading on the top of the gel. To load the sample, use the P20 micropipette with tip to withdraw protein sample, and carefully place the tip on the top of the well. By injecting slowly the sample will sink to the bottom of the well by gravity (**Figure 1**).

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2. If you need to practice, use the P20 micropipette with tip to withdraw 10  $\mu$ L of loading dye from microcentrifuge tube L (pink tube with orange label) on rack. Load the dye into wells 1 to 3 or 7 to 10.

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3. Each of the microcentrifuge tubes M, U1 and U2 (pink tube with orange label) contains 15  $\mu\text{L}$  of protein molecular weight marker, unknown protein U1 and unknown protein U2, respectively. Use micropipette P20 to withdraw 10  $\mu\text{L}$  solution from each tube and load the samples into wells 4 to 6 as shown in **Figure 1**.

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4. As soon as you finish sample loading, **Lift the sign**, lab assistants will connect the power cord to power supply and set the voltage to 200 V for you. The gel will run for 25 minutes. The timer will be set up by an assistant to countdown.

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5. After finishing electrophoresis, **Lift the sign**, lab assistants will disassemble the electrophoresis set-up and give back your gel. Wipe clean the surface of gel with tissue papers and **label the gel with your Student Code** sticker. Lab assistants will take the photo of your gel. Put the photo on the answer sheet using double-sticker (5 points).

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**Answer the following questions:** 回答下列問題

**Q.1.1. (2 points)**

**Figure 2** shows a photograph of a SDS-PAGE gel. The electrophoresis start point and dye-front are indicated. Which side of the gel should be connected to the anode (+ charge) of the power supply? Mark your answer (X) on the answer sheet.

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**Q.1.2. (8 points)**

Based on the information provided in **Figure 2**, make a plot of log molecular weight values of the five marker proteins versus their relative migration- $R_f$  values on the answer sheet (4 points). Use the graph to estimate the molecular weights of unknown proteins on lanes A and B (4 points). Write down your answers on the answer sheet.

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**Q.1.3. (5 points)**

A protein complex of molecular weight 246 kDa is composed of multiple subunits bound by non-covalent interaction. Two protein bands of 57 and 33 kDa were identified after SDS-PAGE. How many 57-kDa and 33-kDa subunits, respectively, are included in the protein complex? Write down your answers on the answer sheet.

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**Q.1.4. (5 points)**

The average molecular weight of amino acid residues is about 110 daltons. How many amino acids are there in the 33-kDa protein subunit? How many nucleotides of RNA are translated into the protein? Write down your answers on the answer sheet.

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**Q.1.5. (5 points)**

Suppose the average molecular weight of nucleotides is 330 daltons. Excluding intron and stop codon, what is the mass ratio of dsDNA encoding the 33-kDa protein to the 33-kDa protein? Write down your answer on the answer sheet.

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**Q.1.6. (5 points)**

Suppose a protein P can bind to a protein Q (MW = 1000 daltons). The binding can be revealed by gel-mobility shift assay. Now 200 pmol of protein P were mixed with various amounts (0 to 500 ng) of protein Q. These mixtures were resolved by 10% (w/v) polyacrylamide gel. Gel was stained by Coomassie blue and is shown in **Figure 3**. Calculate the binding molar ratio of proteins P and Q? Write down your answer on the answer sheet.

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### Protein quantification 蛋白質定量

#### Introduction: 簡介

Coomassie Brilliant Blue G-250 (CBG) is a protein staining reagent. It appears in a different color under different pH conditions. It looks reddish brown in acidic solution, whereas it turns blue under neutral or alkaline condition. Since proteins can provide a relative neutral environment, CBG will turn blue with the maximum absorbance at a wavelength of 595 nm when binding to protein. The more protein there is in a sample, the more CBG will bind to it, and thus, the higher intensity of blue color will be. In other words, the absorbance at 595 nm is proportional to the amount of protein in a sample. Based on this, one can determine the concentration of protein by measuring the blue intensity of a sample.

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#### In the problem set, you will perform the following experiment:

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1. To make BSA concentration standards (**Table 1**), add 0, 2, 4, 6, 8 and 10  $\mu\text{L}$  of 0.5 mg/mL BSA (green color) in A1 to A6 wells of a microplate (**Figure 4**). Make duplicated BSA concentration standards in B1 to B6 wells. If this step is incorrect, you can repeat the procedure in wells A7 to A12 and/or B7 to B12. Adjust the total volume of each BSA solution to 10  $\mu\text{L}$  by adding an appropriate volume of  $\text{H}_2\text{O}$  (**Table 1**).

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2. Add 200  $\mu\text{L}$  of CBG reagent per well in A1 to A6 & B1 to B6. Mix and observe the color change.

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3. To determine the two concentrations X and Y of enzyme E, add various amounts (2, 4, 6, 8 and 10  $\mu\text{L}$ ) of enzyme E (green color) in duplicate to empty wells and bring up the volume to 10  $\mu\text{L}$  with  $\text{H}_2\text{O}$ .

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4. Add 200  $\mu\text{L}$  of CBG reagent per well to the diluted enzyme E. Mix and observe the color change.

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5. **Lift the sign**, lab assistants will accompany you to measure the absorbance values of your samples at 595 nm using spectrophotometer. Put your Student Code on the print-out data with marker pen.

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6. Return to your work bench, and put the result on the answer sheet using double-sticker.

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**Answer the following questions: 回答下列問題**

**Q.2.1. (10 points)**

Calculate the concentrations of BSA in each sample (10  $\mu\text{L}$ ) and fill in the blanks in the table on the answer sheet (Q.2.1.1. 5points). Use these values to plot a standard curve of BSA concentrations (X-axis) versus mean absorbance values of duplicated standards (Y-axis) on the answer sheet (Q.2.1.2. 5points).

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**Q.2.2. (12 points)**

Choose the best sample solution of diluted solution X and Y within the range of BSA standard curve and fill in the table on the answer sheet.

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**Q.2.3. (8 points)** Based on the best sample solution you chose, calculate the original concentrations (X and Y) of enzyme E from the standard curve of BSA concentration. The concentrations should be expressed in units of mg/mL. Write down your answers on the answer sheet.

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### Protein purification

#### **Introduction:**

Column chromatography is commonly used for purification of proteins. The column is made by packing solid porous material (stationary phase) in a column filled with buffer solution (mobile phase). The protein solution to be separated is loaded on top of the column and allowed to percolate into the solid matrix (stationary phase). A reservoir at the top supplies elution buffer constantly which flows through the matrix and passes out of the column at the bottom (the eluent). Since proteins interact with solid matrix in different degree, individual proteins migrate faster or more slowly through the column depending on their properties. Therefore, one can obtain purified proteins by collecting eluent at different times (**Figure 5**).

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Ion-exchange chromatography can be used to separate proteins with different electric charge at a given pH. In anion exchange chromatography, negatively charged proteins bind to positively charged stationary phase. Using solution containing anions to compete with proteins for the adsorption of solid matrix, the bound proteins will be eluted. In practical, proteins are eluted first with buffer containing lower concentration of anion, then with buffer containing higher concentration of anion. Since different charged proteins interact with the stationary phase in different strength, they can be separately eluted by different concentrations of anionic buffers.

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**In the problem set, you will perform the following experiment (5 points):**

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1. Label six 15-mL centrifuge tubes (yellow cap) a1 to a3 and b1 to b3 accordingly, with a marker pen.

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2. Take the anion chromatography column (**Figure 6A**), un-plug the tube and allow the solution to be drained by gravity in the same centrifuge tube. Plug the tube intermediately when the liquid surface reaches the top of the disc (**Figure 6A**, white arrow). Do not over-dry the gel as it may affect protein purification.

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3. Withdraw 200  $\mu$ L of protein solution from microcentrifuge tube C (blue tube with blue label) using a P200 micropipette, apply the sample to the chromatography column slowly by touching the filled pipette tip lightly against the inside wall of the tube (**Figure 6B**).

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4. Un-plug the column and allow the protein sample to drain out, then transfer the column to centrifuge tube a1 (yellow cap). Withdraw 3 mL of anion buffer A (green cap) with a plastic dropper and apply the solution to gel by touching pipette tip against the wall of the tube (**Figure 6C**).

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5. Collect ~1 mL eluent in centrifuge tubes a1 to a3 (yellow cap) sequentially. It takes about 2 to 3 minutes for each tube.

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6. Allow the contents of the column to drain entirely out then transfer the column to centrifuge tube b1 (yellow cap). Withdraw 3 mL of anion buffer B (green cap) with a plastic dropper and apply the solution to gel by touching pipette tip against the wall of the tube (**Figure 6C**).

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7. Collect ~1 mL eluent in centrifuge tubes b1 to b3 (yellow cap) sequentially. It takes about 2 to 3 minutes for each tube.

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
8. Withdraw 50  $\mu$ L of eluent from tubes a1 to a3 & b1 to b3 (yellow cap) and transfer to centrifuge tubes A1 to A3 & B1 to B3 (red cap), respectively. Mix and observe color change. CBG (see introduction in Task II) reagent in tubes A1 to A3 & B1 to B3 will turn blue when it reacts with the eluted protein.

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9. After finishing all the experiments, **Lift the sign**, lab assistants will take photo of your experiment results and put a stamp mark on your answer sheet. Without the stamp mark, you will not be evaluated the Q.3.1.1. and Q.3.1.2.

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**Q.3.1. (7 points)**

Mark the deepest color change (X) on the answer sheet (Q.3.1.1. 5 points). Which of the following buffers (buffer A or buffer B) can be used to elute the protein? Mark your answer (X) on the answer sheet (Q.3.1.2. 2 points).

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**Q.3.2. (5 points)**

Enzyme A is a protein whose surface is evenly distributed with electric charges. If enzyme A can be eluted from anionic exchange chromatography by high concentration of anionic buffer, what is the property of enzyme A with respect to electric charge? Mark (X) the answer on the answer sheet.

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
**Q.3.3. (4 points)**

Different amino acids differ in the chemical nature of the R group (side chain). **Figure 7** shows four amino acids A, B, C, and D in their prevailing ionic forms at pH 7.2, with the side chain marked in white box. Which of the following amino acids in **Figure 7** would be present more frequently on enzyme A? Write down your answer (X) on the answer sheet.

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**Figure 7**

**Q.3.4. (5 points)**

Hydrophobic interaction chromatography can be used to separate proteins based on their hydrophobicity. To perform the chromatography, protein samples were first treated with buffer containing high concentration of salts such as ammonium sulfate  $(\text{NH}_4)_2\text{SO}_4$ , which will remove water molecules from the protein surface. This causes the hydrophobic area on the surface of the protein to be exposed. When the salt-treated proteins are subjected for chromatography, they will be absorbed on the stationary phase through hydrophobic interactions. The higher the hydrophobicity of the protein, the stronger the absorption. As salt concentration can affect the hydrophobic interaction between the protein and the stationary phase, different proteins can be separately eluted by using different concentrations of salt-containing buffers. If enzyme A is highly hydrophobic, which of the following buffers should be used to separate enzyme A from other proteins by chromatography? Mark (X) the answer on the answer sheet.

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(E) High-salt buffer first then low-salt buffer

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**Q.3.5. (4 points)**

If enzyme A is highly hydrophobic, which of the amino acids in **Figure 7** would be present more frequently on enzyme A? Mark (X) the answer on the answer sheet.

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**Q.3.6. (5 points)**

Gel filtration chromatography separates proteins based on their sizes. The gel, or stationary phase, consists of cross-linked polymer beads with engineered pores of a particular size. Small proteins enter the pores and are retarded by their more labyrinthine path. Large proteins cannot enter the pores and so take a short path through the column, around the beads. **Table 2** is a list of gels and their fractionation ranges. Suppose both enzyme A (22 kDa) and protein B (44 kDa) are single-subunit proteins. One would like to purify enzyme A from a mixture containing enzyme A and protein B using gel filtration chromatography. Which gel is best suited for the job? Mark your answer (X) on the answer sheet.

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**Q.3.7. (5 points)**

Assume that the concentration of total proteins in the original solution is 1 mg/mL and the activity of enzyme A is 0.5 units in 1-mL protein sample. The concentration of total proteins after purification is 0.1 mg/mL and the activity of enzyme A is 1 unit in 1-mL protein sample. Calculate the purification factor (times of purity improvement) of enzyme A. Write down your answer on the answer sheet.

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