Country: _____

Student Code: _____

23rd INTERNATIONAL BIOLOGY OLYMPIAD

8th – 15th July, 2012

SINGAPORE

PRACTICAL TEST 2

MICROBIOLOGY & BIOCHEMISTRY

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Total points: **100** 總分：100

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Dear Participants

- In this test, you have been given the following two tasks:

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Task 1: Bacteriophage: an effective agent in the killing of bacteria. (50 points)

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Part A: Effects of Phage and antibiotics on the killing of antibiotic-resistant *E. coli* (31 points)

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Task 2: Titration of an amino acid. (50 points)

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- Use the Answer Sheet, which is provided separately, to answer all the questions.

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- The answers written in the Question Paper will **NOT** be evaluated.

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- Please make sure that you have received all the materials and equipment listed for each task.

If any of these items are missing, please raise your hand **immediately**.

- Stop answering and put down your pen IMMEDIATELY when the bell rings.

- At the end of the test, place the Answer Sheets and Question paper in the envelope provided.

Our Assistants will collect the envelope from you.

Have fun and Good Luck! 😊

Materials and equipment: 材料和器材

For Task I: Bacteriophage: an effective agent in the killing of bacteria

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Task I (50 points)

Bacteriophage: an effective agent in the killing of bacteria

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Part A. Effects of Phage and antibiotics on the killing of antibiotic-resistant *E. coli* (31 points)

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Introduction

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A bacteriophage is a virus that infects bacterial cells. Certain bacteriophages can kill bacteria cells by lysis. The bacteriophage is now recognized as an effective agent in the killing of pathogenic bacteria. This provides a good alternative to antibiotics in our combat against disease-causing bacteria that might be resistant to traditional antibiotics.

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You are required to design a simple experiment, with proper controls, to examine the killing efficiency of phage of an ampicillin-resistant *E. coli*. Answer the following questions **in the Answer Sheet** and follow the instructions given below.

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Q1.1 (1 point) To dilute the *E. coli* culture from 1×10^7 cells/ml to 2×10^5 cells/ml, what would be the dilution factor needed?

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Q1.2 (1 point) For 1 ml of *E coli* culture at a cell density of 2×10^5 cells/ml, the final concentration of ampicillin used should be 10 μ g/ml. What would be the volume of ampicillin stock (1mg/ml) used?

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Q1.3 (1 point) For 1 ml of *E coli* culture at a cell density of 2×10^5 cells/ml, final titre of phage used should be. What would be the volume of phage stock (10^8 pfu/ml) used?

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Q1.4 (1 point \times 15 = 15 points) With the above calculated dilution factors, fill in the table **in the Answer Sheet** with your experimental plan. One example is already given for Tube 1 in the table provided. All units are in μ l. Carry out your proposed experiment by incubating the four tubes (placed in the labelled floating rack) for 40 minutes (stop watch provided) in the 37 $^{\circ}$ C water bath assigned to you. Hand over the floating rack to the technician at the water bath.

After incubation, transfer your samples to the cuvettes labelled 1 to 4. In order to observe the killing of bacteria cells, measure the absorbance at 595 nm wavelength. Bring your samples to the spectrophotometer that is allocated for your use and hand over your samples to the technician. You are to record your own readings as the samples are measured.

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Q1.5 (0.75 × 2 + 1.5 points × 6 = 10.5 points) Fill in the absorbance reading at 595 nm of the different tubes of reactions in the table provided **in the Answer Sheet**. Taking 1 absorbance unit of the *E. coli* cells at 595 nm to be equivalent to 1×10^7 cells/ml, what are the cell densities of the *E. coli* in the respective reaction tubes?

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Q1.6 (0.5 points × 5 = 2.5 points) Which of the following are correct? Indicate correct answer(s) with a tick(✓) and incorrect answer(s) with a cross (✗).

- Due to the ampicillin resistance, the bacteria cell wall prevented easy penetration of the antibiotics, but allowed the phages to enter the *E. coli* cells to cause lysis.
- The ampicillin resistance in the *E. coli* did not prevent the ability of the phage to adsorb onto the bacteria cells.
- The bacteriophage likely has a lytic life cycle of around 20 to 30 mins and hence lysis of the *E. coli* was observable during the short experiment.
- The temperature of 37 °C was not the correct temperature for ampicillin to kill the *E. coli* cells.
- The phages competed with the *E. coli* for the nutrients in the LB broth and the bacteria cells lysed due to insufficient nutrients.

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The Table below shows the legends for photographs of *E. coli* lawns that are untreated and infected with bacteriophages. The *E. coli* culture used had a starting cell density of 0.5×10^4 cells/ml. 0.5 ml of phage was used to infect the *E. coli* cells. Serial dilutions of the phage culture were made as indicated and used for infection. (The photographs labelled A to H will be provided as part of materials for the lab task).

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Q1.7 (2 points × 4 = 8 points) Based on the number of plaques observed in the photos, calculate the number of plaques that would be observed if the original undiluted phage culture were used.

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Q1.8 (3 points) To estimate the titre of a phage culture, serial dilutions as shown in the photos (A to H) are normally performed. Based on the number of plaques shown, indicate with a tick (✓) which is the best dilution to confirm the phage titre.

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Q1.9 (4 points × 2 = 8 points) Using the information given and your answers above, determine:

- the plaque forming units per milliliter (pfu/ml) of the phage culture used and
- the multiplicity of infection (defined as the ratio of phages to *E. coli*) at the best dilution determined in Q1.8..

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Titration of an Amino Acid

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Introduction 背景知識

Amino acids are organic molecules possessing both carboxyl and amino groups. Table 1 shows the 20 amino acids that cells use to build their thousands of proteins.

The majority of the standard amino acids are diprotic molecules since they have two dissociable protons: one on the amino group and the other on the carboxyl group; there is no dissociable proton in the R group.

Recall: For an acid HA, the acid dissociation constant for the equilibrium of $HA \rightleftharpoons H^+ + A^-$ is K_a .

$$K_a = [H^+][A^-] / [HA]$$

More often, the strength of acids is expressed in terms of the pK_a of the acid: $pK_a = -\log K_a$

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In the titration of such a diprotic amino acid, the titration will thus occur in two steps as the more acidic carboxyl group (lower pK_{a1}) and the less acidic amino group (pK_{a2}) successively lose their protons.

In addition, the pH at which the net charge on the molecule is zero is called the isoelectric point (pI) of the molecule, a useful constant in characterizing and purifying molecules. Using a titration curve,

the pI can be empirically determined as the inflection point between the pK_a's of the anionic and cationic forms.

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The apparent pK_a values for the two dissociation steps may be extrapolated from the midpoints of each step. This can be shown by the Henderson-Hasselbach equation:

$$\text{pH} = \text{pK}_a + \log \left\{ \frac{[\text{A}^-]}{[\text{HA}]} \right\}$$

The pK_{a1} (pK_a for the carboxyl acid group) is where half the acid group has been titrated. Therefore the equation becomes: $\text{pH} = \text{pK}_a$

Similarly, the pK_{a2} (pK_a for the amino group) can be determined.

In this experiment, you will titrate an unknown amino acid Z and determine its pI, pK_{a1} and pK_{a2}.

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Procedure 實驗步驟

1. Fill the burette with the standardized NaOH solution. Record the exact concentration of this standardized NaOH solution **in the Answer Sheet**.

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2. Pipette 25 ml of the unknown amino acid solution Z into a clean 100 ml beaker.

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3. Carefully place the pH probe and a magnetic stirring bar into the amino acid solution, so that the probe is far enough into the solution, but not touching the stirring bar or beaker. Clamp and adjust the pH probe such that the stirring bar will not hit the probe while stirring. **DO NOT TOUCH THE CALIBRATION.**

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4. Titration 1 滴定 1

Rinse the pH probe with deionized water. Dry the probe gently with a piece of Kimwipe paper.

Determine the pH of the amino acid solution Z before the addition of NaOH.

Next, titrate the amino acid solution with the NaOH from the burette. Add approximately 1.00 ml of the NaOH to the amino acid at a time. Record the exact volume dispensed and the pH of the solution after every 1.00 ml interval **in the Answer Sheet.**

Continue until approximately 25 ml of NaOH has been added.

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Rinse the pH probe with deionized water.

Dry the probe gently with a piece of Kimwipe paper.

Refill the burette with the standardized NaOH solution and repeat steps 2 – 4.

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Q2.1 (3 points × 3 = 9 points) Table 1 shows the chemical structures of the twenty standard amino acids. With reference to these structures, draw structures to show the complete dissociation of glycine, proline and asparagine.

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Q2.2 (3 points × 2 = 6 points) For both titrations, record the volume of NaOH (ml) added during the titration and the observed pH value for the unknown amino acid.

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Q2.3 (5 points × 2 = 10 points) Using your data, plot the graphs of each titration run (pH versus Vol. of NaOH (ml)) in Graphs 1 and 2 provided in the Answer Sheet.

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Q2.4 (2 points × 2 = 4 points) From your titration curves, find the pI and label it on each graph.

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Q2.5 (4 points × 2 = 8 points) Find and label the pK_{a1} and pK_{a2} on each graph.

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Q2.5.1 (2 points × 2 = 4 points) What is the mean pK_{a1} and pK_{a2} ?

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Q2.6 (5 points) 0.9210 g of the unknown amino acid Z was dissolved in 80 ml of deionized water.

Determine the molecular weight of the unknown amino acid Z.

Note: In order to start with a fully protonated amino acid, HCl solution has been added. This is equivalent to 3.2 ml of the NaOH solution. To determine the actual number of moles of NaOH needed to reach the pI, subtract 3.2 ml from the volume of NaOH used to reach the first end point.

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Q2.7 (2 points) Based on Table 2, identify amino acid Z.

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Table 1. Structures of the 20 standard amino acids.

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