25th INTERNATIONAL BIOLOGY OLYMPIAD

5 – 13 July, 2014

INDONESIA

PRACTICAL TEST 1 實作測驗 1

CELL & MOLECULAR BIOLOGY

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Total points: 64.5 總分：64.5

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
COUNTRY:
STUDENT:

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- In this test, you have been given the following tasks:
本測驗中，你有以下實驗要完成
 - Part A.** Confirmation of plasmid samples X, Y, and Z by restriction enzyme analysis.
(40 points)
以限制酶切位分析去確認質體 X, Y 和 Z
 - Part B.** Cell reproduction and telomere analysis of *Paramecium*. (24.5 points)
草履蟲細胞生殖和端粒分析
- Answer all the questions in the **Answer Sheet** provided.
所有答案要寫在答案卷上
- The answers written in the **Question Paper** will **NOT** be evaluated.
寫在試題卷上的答案將不給分
- Write your answers legibly in ink.
使用原子筆清楚寫答案
- Please make sure that you have received all the materials and equipment listed for each task. **If any of these items are missing, please raise your hand immediately.**
確認所有的材料、儀器都符合清單所列，若發現任何項目缺少，立即舉手
- **Use the gloves provided for the experiment**
實驗時使用提供的手套
- Stop answering and put down your pen **immediately** when the final bell rings.
測驗結束鈴響時，立刻放下筆停止作答
- At the end of the test, place the Answer Sheet and Question Paper in the envelope provided. Our Assistants will collect the envelope from you.
測驗結束時，將答案卷和試題卷放進提供的信封套中，助理會來收取信封套
- **Note!** All the DNA electrophoresis instruments will be turned on simultaneously by the assistants **50 minutes after the start of the test**. Make sure that you have placed your samples in the agarose gel in accordance to the instructions on the Question paper before the power supply is turned on. After this, you will **NOT** be allowed to run the gel. To remind you, the assistant will ring the bell three times. Bell 1 will remind you that the electrophoresis will be carried out in 20 minutes. Bell 2 is a warning that electrophoresis will begin in two

minutes' time. Bell 3 will mark the end of electrophoresis. At that time you will be asked to remove your gel and place it into the box provided for assessment.

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Note : Use given reagents properly! No additional reagents will be provided.

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Part A (40 points)

Identification of Plasmids by Restriction Enzyme Analysis

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Introduction 背景介紹

A scientist was surprised to find unlabelled tubes in the freezer. These tubes, containing Plasmids X, Y and Z, were arbitrarily labelled Plasmid 1, Plasmid 2 and Plasmid 3. These plasmids are indistinguishable by DNA electrophoresis because they are all 3750 bp long. Since the three plasmids differ in their restriction pattern when using *EcoRI* and/or *HindIII* (Figure 1), you will conduct restriction enzyme analysis to correctly identify the three plasmids.

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
Figure 1. Restriction maps of Plasmids X, Y, Z.

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Table I in the **Answer Sheet** shows the experimental design you will perform. The experiment is divided into Series 1 and Series 2, carried out simultaneously. In each series, all three plasmids should be cut with the same restriction enzyme(s). The final concentration of the restriction buffer should be 1X.

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Question 1.1. Determine how these restriction enzyme(s) can be used to differentiate between the plasmids. Complete Table I on the **Answer Sheet to indicate your chosen protocol for each series!** (8 points). You must choose one or two enzymes for each series that are most informative based on Figure 1. If the enzyme is used in the reaction, add 1 μL enzyme.

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Question 1.2. Complete the figure on your **Answer Sheet** to show the expected molecular weight(s) of the completely digested plasmid for Plasmids X, Y and Z that you would observe using each of your two chosen Series. Indicate the restriction enzyme(s) used for Series 1 and Series 2 of your experimental design. (6 points – 1 points each lane).

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Plasmid restriction and DNA electrophoresis protocol

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Note: Electrophoresis will start 50 minutes after the test begins. After this, you will NOT be allowed to run the gel.

Note: Spin down all reagents in microtubes before pipetting (Be sure to balance the minicentrifuge by placing the microtubes opposite each other, even if there is only one tube to spin).

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1. Label the microtubes S1 to S6, and prepare the reaction mixtures according to Table I in the Answer Sheet. Spin down the mixture by placing the microtubes in the minicentrifuge for 10-20 seconds.

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2. Place the microtubes into the floating rack labelled with your bench number and incubate for at least 10 minutes at 37°C in the water bath, located at the end of your aisle in the direction of the arrow.

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3. Add **ONLY 2 μ L** of DNA staining solution into microtube labelled 1 Kb Ladder. Spin down the solution for 10 seconds.

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4. After incubation, raise your bench number and your tubes will be returned to you. Add **ONLY** **1 μ L** of DNA staining solution into each of these tubes, mix well by pipetting and spin down any residual liquid.

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5. Open the gel agarose package and place the gel with its tray in the electrophoresis chamber. Pour the entire TAE buffer from the bottle into the tank.

Note: Black plug is cathode (-) and red plug is anode (+).

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6. Load 10 μ l of samples S1 to S6 and molecular weight marker solution into the agarose gel as shown in **Figure 2**.

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Figure 2 圖 2

7. Close the lid of the gel chamber and notify the assistant by raising your hand that you are ready. The assistant will connect the cables to the power supply.

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During the electrophoresis continue with Part B.

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8. After the electrophoresis (indicated by Bell 3 ringing), the assistants will turn off the power supply. Carefully place the gel with the tray in the plastic box provided. **Pour the running buffer from the electrophoresis chamber into the plastic box.** Close the box, label the sticker with your **Student ID** and affix it to the side of the box. Leave the box on your bench. Later the gel will be documented and the picture will be attached onto your answer sheet by the assistants.

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Part B (24.5 points)

Cell Reproduction and Telomere Analysis of *Paramecium*

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Two of the phenomena displayed by *Paramecium tetraurelia* are binary fission (asexual reproduction) and conjugation (exchange of genetic materials), which are affected by the abundance of nutrients in the media. Figure 3 shows microscopic observations of two *P. tetraurelia* cultures grown in rich and poor media, respectively.

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Question 2.1. (2 points) Which of the phenomena are shown by the two cultures in Figure 3A and 3B?

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
Figure 3. *P. tetraurelia* growing in rich (A) and poor (B) media.

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Paramecium, which has just undergone conjugation, continues to asexually. In this instance, telomere shortening occurs for a certain number of divisions during logarithmic growth. **Table B-1** shows data for the growth of a strain of *Paramecium* in three replicates for four days. The cultures can be grown continuously with no significant change in the fission rate.

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Table B-1. *P. tetraurelia* growth in the first four days.

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Question 2.2. (2.5 + 8 points) Calculate the average cell concentrations and using the logarithmic value, draw the growth curve of *P. tetraurelia* from Days 0 to 4 on the graph provided in your **Answer Sheet** (Precision: Integer).

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Figure 4 shows the result of Southern blot analysis of *P. tetraurelia* telomeres for a total of 30 synchronized generations of cultivation. DNA from *P. tetraurelia* was digested at the beginning of the telomere region, blotted onto a membrane and then hybridized with telomere probes. Telomeres were seen as smears, the midpoints of which can be matched with molecular weight markers to determine the average telomere length.

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
Figure 4. Southern blot of *P. tetraurelia* telomeres. Lanes 1 to 5 show telomeres after 4, 7, 17, 23, and 30 generations, respectively. DNA size markers (M) are shown next to Lane 5 (in kbp)

Internal: Internal part of the chromosomal DNA that contains telemetric sequences)

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Question 2.3. (12 points) Indicate in the answer sheet with a tick (✓) if the following statements are true or false.

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Question 2.3. (12 points)

No.	Statement 敘述
a	The media can support the growth of <i>P. tetraurelia</i> to more than 10^4 cells/mL. 培養基能提供草履蟲生長到大於 10^4 cells/mL
b	The Telomere probe used for Southern blot analysis is not specific for telomeric sequences. 用於南方墨點法的端粒探針，對於端粒序列的專一性不足
c	Based on the Southern blot, telomeres of the cells from the same generation have the same length. 根據南方墨點法的結果，同一世代的端粒具有相同的長度
d	Paramecium can grow from 10^5 to 10^7 cells in less than three days, if enough medium is provided. 如果營養源供應充足，草履蟲能在三天內從 10^5 長到 10^7
e	Based on telomere data from lane 4 and 5, telomeres shorten by 20-25 base pairs per generation. 根據第 4 與 5 行的端粒資料，每個世代的端粒會縮短約 20-25 base pairs
f	At beginning of the experiment (parental generation), telomeres were likely between 1500 and 1700 base pairs long. 在實驗開始進行前（親代），端粒長度約在 1500 - 1700 base pairs

25th INTERNATIONAL BIOLOGY OLYMPIAD

5 – 13 July, 2014

INDONESIA

PRACTICAL TEST 1 實作測驗 1

CELL & MOLECULAR BIOLOGY 細胞與分子生物學

ANSWER SHEET 答案卷

Total points: **64.5** 總分： **64.5**

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COUNTRY:
STUDENT ID:

Task (64.5 points)

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Plasmid Identification and Telomere Analysis

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- The answers have to be given either with a tick (✓) or with Arabic numbers. The numbers "1" and "7" can look very similar in handwriting. To make sure that those two numbers can be well distinguished by the IBO staff, please write them as you normally would into the following box.
答案必須以勾選 (✓) 或 數字 的方式回答，手寫數字 "1" 和 "7" 常被混淆，為確保試務人員正確判讀你的答案，請在下列方框中寫下你習慣的 "1" 和 "7" 手寫樣式

1 =		7 =	
------------	--	------------	--

Part A. Confirmation of plasmid sample X, Y, and Z by restriction analysis and DNA electrophoresis. (40 points)

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Q 1.1. (8 points)

Table I. Design of Experiment for Plasmid Identification

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* If the enzyme is used in the reaction, add 1 μL enzyme

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
Q 1.2. (6 points)

Please give a tick (✓) to the chosen enzyme on the table below and put minus (-) if the enzyme is not chosen.

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
Enzyme Used	Series 1	Series 2
<i>EcoR1</i>		
<i>HindIII</i>		

Q 1.3. Photo (26 points)

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Part B. Cell reproduction and telomere analysis of *Paramecium* (24.5 points)

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Q 2.1. (2 points).

Please give a tick (✓) for the presence and minus (-) for the absence of the phenomenon

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